

Product Datasheet

Mouse TGFB3 ELISA Kit (orb409559)

Description	Mouse Transforming Growth factor β3,TGF-β3 ELISA kit
Reactivity	Mouse
Range	125 pg/mL-8000 pg/mL
Target	TGFB3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P17125
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.