

Product Datasheet

Pig AFP ELISA Kit (orb409577)

Description	Pig Alpha-fetoprotein(AFP) ELISA kit
Reactivity	Porcine
Range	3.12 ng/mL-200 ng/mL
Target	AFP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q8MJ76
Sensitivity	0.78 ng/mL
Expiration Date	Please enquire.