

Product Datasheet

Pig CD163 ELISA Kit (orb409652)

Description	Pig Scavenger receptor cysteine-rich type 1 protein M130(CD163) ELISA kit
Reactivity	Porcine
Range	93.75 ng/mL-6000 ng/mL
Target	CD163
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q2VL90
Sensitivity	23.4 ng/mL
Expiration Date	Please enquire.