

Product Datasheet

Pig IFN β /IFNB ELISA Kit (orb409686)

Description	Pig Interferon β , IFN- β /IFNB ELISA Kit
Reactivity	Porcine
Range	50 pg/mL-2000 pg/mL
Target	IFN- β /IFNB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q29412
Sensitivity	25 pg/mL
Expiration Date	Please enquire.