

Product Datasheet

Pig IL15 ELISA Kit (orb409694)

Description	Pig Interleukin 15(IL-15) ELISA Kit
Reactivity	Porcine
Range	0.78 pg/mL-50 pg/mL
Target	IL15
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q95253
Sensitivity	0.195 pg/mL
Expiration Date	Please enquire.