

Product Datasheet

Pig IL1 alpha ELISA Kit (orb409697)

Description	Pig Interleukin 1 α ,IL-1 α ELISA kit
Reactivity	Porcine
Range	25 pg/mL-1600 pg/mL
Target	IL1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P18430
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.