


Product Datasheet

CYP3A43 antibody (orb40972)

Description	CYP3A43 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cytochrome P450 3A43 protein (244-503AA)
Target	CYP3A43
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HB55
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver(lan...