

Product Datasheet

Pig MMP-9 ELISA Kit (orb409753)

Description	Pig matrix metalloproteinase 9,MMP-9 ELISA Kit
Reactivity	Porcine
Range	0.625 ng/mL-40 ng/mL
Target	MMP-9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q2VI03
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.