

Product Datasheet

Pig TIMP1 ELISA Kit (orb409754)

Description	Pig tissue inhibitors of metalloproteinase 1, TIMP-1 ELISA Kit
Reactivity	Porcine
Range	7.8 ng/mL-500 ng/mL
Target	TIMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35624
Sensitivity	1.95 ng/mL
Expiration Date	Please enquire.