

## Product Datasheet

### Pig BNP ELISA Kit (orb409776)

|                          |                                                 |
|--------------------------|-------------------------------------------------|
| <b>Description</b>       | Pig brain natriuretic peptide, BNP<br>ELISA Kit |
| <b>Reactivity</b>        | Porcine                                         |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                          |
| <b>Target</b>            | BNP                                             |
| <b>Note</b>              | For research use only                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates               |
| <b>Uniprot ID</b>        | <b>P07634</b>                                   |
| <b>Sensitivity</b>       | 7.8 pg/mL                                       |
| <b>Expiration Date</b>   | Please enquire.                                 |