

Product Datasheet

Pig NOS3 ELISA Kit (orb409778)

Description	Pig Endothelial nitric oxide synthase,eNOS ELISA Kit
Reactivity	Porcine
Range	3.12 IU/mL-200 IU/mL
Target	NOS3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q28969
Sensitivity	0.78 IU/mL
Expiration Date	Please enquire.