

Product Datasheet

INPP4A antibody (orb40981)

Description	INPP4A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Type I inositol 3,4-bisphosphate 4-phosphatase protein (1-250AA)
Target	INPP4A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96PE3
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.