

Product Datasheet

Pig ANGPT2 ELISA Kit (orb409833)

Description	Pig Angiopoietin 2,ANG-2 ELISA Kit
Reactivity	Porcine
Range	0.94 ng/mL-60 ng/mL
Target	ANGPT2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BDY7
Sensitivity	0.235 ng/mL
Expiration Date	Please enquire.