

Product Datasheet

Pig Chlamydia suis antibody (IgG) ELISA Kit [Out of stock] (orb409842)

Description	Pig chlamydia suis antibody (IgG) ELISA kit
Reactivity	Porcine
Range	Request Information
Tested Applications	ELISA
Target	C. suis Ab (IgG)
Note	For research use only
Application notes	Sample Type: serum, plasma Assay Time: 1-3 hours Sample Volume: 50-100 µl Detection Wavelength: 450nm
Sample Types	serum, plasma
Sensitivity	Request Information
Expiration Date	Please enquire.