

Product Datasheet

Pig NGAL ELISA Kit (orb409864)

Description	Pig neutrophil gelatinase-associated lipocalin(NGAL)ELISA Kit
Reactivity	Porcine
Range	7.8 pg/mL-500 pg/mL
Target	NGAL
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	F1RRX1
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.