

## Product Datasheet

### Rabbit ORM1 ELISA Kit (orb409877)

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Description</b>       | Rabbit Alpha-1-acid glycoprotein<br>1(ORM1) ELISA kit |
| <b>Reactivity</b>        | Rabbit                                                |
| <b>Range</b>             | 15.6 ng/ml-1000 ng/ml                                 |
| <b>Target</b>            | ORM1                                                  |
| <b>Note</b>              | For research use only                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                     |
| <b>Uniprot ID</b>        | <b>P25227</b>                                         |
| <b>Sensitivity</b>       | 3.9 ng/ml                                             |
| <b>Expiration Date</b>   | Please enquire.                                       |