

Product Datasheet

MOCS1 antibody (orb40992)

Description	MOCS1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Molybdenum cofactor biosynthesis protein 1 protein (90-360AA)
Target	MOCS1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZB8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart tis...