

Product Datasheet

Rabbit SFTPB ELISA Kit (orb409921)

Description	Rabbit Pulmonary surfactant-associated protein B(SFTPB) ELISA kit
Reactivity	Rabbit
Range	1.87 ng/mL-120 ng/mL
Target	SFTPB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P15285
Sensitivity	0.46 ng/mL
Expiration Date	Please enquire.