

Product Datasheet

Rabbit MMP2 ELISA Kit (orb409930)

Description	Rabbit matrix metalloproteinase 2/Gelatinase A,MMP-2 ELISA Kit
Reactivity	Rabbit
Range	1.56 ng/mL-100 ng/mL
Target	MMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P50757
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.