

Product Datasheet

Rabbit SDF-1/CXCL12 ELISA Kit (orb409931)

Description	Rabbit stromal cell derived factor 1 (SDF-1/CXCL12) ELISA kit
Reactivity	Rabbit
Range	0.57 ng/mL-10 ng/mL
Target	SDF-1/CXCL12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Sensitivity	0.36 pg/mL
Expiration Date	Please enquire.