

Product Datasheet

Rabbit Apolipoprotein AI ELISA Kit (orb409962)

Description	Rabbit apolipoprotein A1 (Apo-A1) ELISA Kit
Reactivity	Rabbit
Range	6.25 µg/mL-400 µg/mL
Target	APOA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09809
Sensitivity	1.56 µg/mL
Expiration Date	Please enquire.