

Product Datasheet

MSRB3 antibody (orb40997)

Description	MSRB3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Methionine-R-sulfoxide reductase B3 protein (1-185AA)
Target	MSRB3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IXL7
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.