

Product Datasheet

Rabbit MCP1 ELISA Kit (orb410003)

Description	Rabbit monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP1/MCAF ELISA kit
Reactivity	Rabbit
Range	23.5 pg/mL-1500 pg/mL
Target	CCL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P28292
Sensitivity	5.87 pg/mL
Expiration Date	Please enquire.