

Product Datasheet

Rabbit BGLAP ELISA Kit (orb410010)

Description	Rabbit Osteocalcin/bone gla protein,OT/BGP ELISA Kit
Reactivity	Rabbit
Range	1.25 ng/mL-80 ng/mL
Target	BGLAP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P39056
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.