

## Product Datasheet

### Rabbit BALP ELISA Kit (orb410011)

**Description** Rabbit Bone Alkaline Phosphatase (BALP) ELISA Kit

**Reactivity** Rabbit

**Range** 0.156 mIU/mL-10 mIU/mL

**Target** BALP

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma

**Uniprot ID** **G1SLU6**

**Sensitivity** 0.039 mIU/mL

**Expiration Date** Please enquire.