

Product Datasheet

Rabbit MMP13 ELISA Kit (orb410018)

Description	Rabbit Matrix Metalloproteinase 13,MMP-13 ELISA Kit
Reactivity	Rabbit
Range	93.75 pg/mL-6000 pg/mL
Target	MMP13
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O62806
Sensitivity	23.44 pg/mL
Expiration Date	Please enquire.