

Product Datasheet

Rabbit TIMP1 ELISA Kit (orb410021)

Description	rabbit tissue inhibitors of metalloproteinase 1,TIMP-1 ELISA Kit
Reactivity	Rabbit
Range	15.6 pg/mL-1000 pg/mL
Target	TIMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates , tissue homogenates
Uniprot ID	P20614
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.