

Product Datasheet

Rabbit TRAP5b ELISA Kit (orb410025)

Description	Rabbit Tartrate-resistant acid phosphatase 5b (TRAP5b) ELISA Kit
Reactivity	Rabbit
Range	78 µU/mL-5000 µU/mL
Target	TRAP5b
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O97860
Sensitivity	19.5 µU/mL
Expiration Date	Please enquire.