

Product Datasheet

rabbit ICAM-1/CD54 ELISA Kit (orb410044)

Description rabbit intercellular adhesion molecule
1,ICAM-1/CD54 ELISA Kit

Reactivity Rabbit

Range 0.625 ng/mL-40 ng/mL

Target ICAM-1/CD54

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma

Uniprot ID **U3KPG6**

Sensitivity 0.156 ng/mL

Expiration Date Please enquire.