

Product Datasheet

Rabbit TGFB1 ELISA Kit (orb410059)

Description	rabbit Transforming Growth factor β 1, TGF- β 1 ELISA Kit
Reactivity	Rabbit
Range	2.86 ng/mL-50 ng/mL
Target	TGFB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q95215
Sensitivity	1.78 ng/mL
Expiration Date	Please enquire.