

## Product Datasheet

### Rat APP ELISA Kit (orb410098)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Description</b>       | Rat Amyloid beta A4 protein(APP)<br>ELISA kit |
| <b>Reactivity</b>        | Rat                                           |
| <b>Range</b>             | 15.6 ng/mL-1000 ng/mL                         |
| <b>Target</b>            | APP                                           |
| <b>Note</b>              | For research use only                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates             |
| <b>Uniprot ID</b>        | <b>P08592</b>                                 |
| <b>Sensitivity</b>       | 3.9 ng/mL                                     |
| <b>Expiration Date</b>   | Please enquire.                               |