

Product Datasheet

Rat ANGPTL4 ELISA Kit (orb410100)

Description	Rat Angiopoietin-related protein 4(ANGPTL4) ELISA kit
Reactivity	Rat
Range	7.8 pg/mL-500 pg/mL
Target	ANGPTL4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, bronchoalveolar lavage fluid, urine
Uniprot ID	Q6TMA8
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.