

Product Datasheet

Rat ARRB2 ELISA Kit (orb410123)

Description	Rat Beta-arrestin-2(ARRB2) ELISA kit
Reactivity	Rat
Range	46.88 pg/mL-3000 pg/mL
Target	ARRB2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P29067
Sensitivity	11.72 pg/mL
Expiration Date	Please enquire.