

## Product Datasheet

### Rat CHI3L1 ELISA Kit (orb410170)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| <b>Description</b>       | Rat Chitinase-3-like protein<br>1(CHI3L1) ELISA kit |
| <b>Reactivity</b>        | Rat                                                 |
| <b>Range</b>             | 4.7 ng/ml - 300 ng/ml                               |
| <b>Target</b>            | CHI3L1                                              |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                        |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                   |
| <b>Uniprot ID</b>        | <b>Q9WTV1</b>                                       |
| <b>Sensitivity</b>       | 1.17 ng/ml                                          |
| <b>Expiration Date</b>   | Please enquire.                                     |