

Product Datasheet

Rat CX3CR1 ELISA Kit (orb410208)

Description	Rat CX3C chemokine receptor 1(CX3CR1) ELISA kit
Reactivity	Rat
Range	23.5 pg/mL-1500 pg/mL
Target	CX3CR1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P35411
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.