

Product Datasheet

Rat FNDC5 ELISA Kit (orb410257)

Description	Rat Fibronectin type III domain-containing protein 5(FNDC5) ELISA kit
Reactivity	Rat
Range	0.625 ng/mL-40 ng/mL
Target	FNDC5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8K3V5
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.