

Product Datasheet

Rat GUCA2B ELISA Kit (orb410310)

Description	Rat guanylate cyclase activator 2B (uroguanylin) (GUCA2B) ELISA kit
Reactivity	Rat
Range	7.8 pg/mL-500 pg/mL
Target	GUCA2B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P70668
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.