

Product Datasheet

Rat Hepcidin-25 ELISA Kit (orb410321)

Description	Rat Hepcidin(HAMP) ELISA kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	HAMP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q99MH3
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.