

Product Datasheet

Rat CSF1R ELISA Kit (orb410379)

Description	Rat Macrophage colony-stimulating factor 1 receptor(CSF1R) ELISA kit
Reactivity	Rat
Range	0.312 ng/mL-20 ng/mL
Target	CSF1R
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q00495
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.