

Product Datasheet

Rat PTGDS ELISA Kit (orb410458)

Description	Rat Prostaglandin-H2 D-isomerase(PTGDS) ELISA kit
Reactivity	Rat
Range	9.37 ng/mL-600 ng/mL
Target	PTGDS
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cerebrospinal fluid (CSF), tissue homogenates
Uniprot ID	P22057
Sensitivity	2.34 ng/mL
Expiration Date	Please enquire.