



Product Datasheet

DCTN6 antibody (orb41047)

Description	DCTN6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Dynactin subunit 6 protein (1-190AA)
Target	DCTN6
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000399
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(lane 1), U...



Immunohistochemical staining of human ga...