

Product Datasheet

Rat TNNI3K ELISA Kit (orb410502)

Description	Rat Serine/threonine-protein kinase TNNI3K(TNNI3K) ELISA kit
Reactivity	Rat
Range	18.75 pg/mL-1200 pg/mL
Target	TNNI3K
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7TQP6
Sensitivity	4.69 pg/mL
Expiration Date	Please enquire.