

Product Datasheet

Rat ATP1A1 ELISA Kit (orb410509)

Description	Rat Sodium/potassium-transporting ATPase subunit alpha-1(ATP1A1) ELISA kit
Reactivity	Rat
Range	31.2 pg/mL-2000 pg/mL
Target	ATP1A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P06685
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.