



Product Datasheet

EIF3M antibody (orb41055)


Description	EIF3M antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit M protein (1-374AA)
Target	EIF3M
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7L2H7
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), AraC...



Immunohistochemical staining of human placental tissue...



Immunoprecipitation analysis of HeLa cell lysates...