

Product Datasheet

Rat PGP9.5 ELISA Kit (orb410551)

Description	Rat Ubiquitin carboxyl-terminal hydrolase isozyme L1(UCHL1) ELISA kit
Reactivity	Rat
Range	0.78 ng/mL-50 ng/mL
Target	UCHL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q00981
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.