

Product Datasheet

Rat IP3 ELISA Kit (orb410567)

Description	Rat inositol 1,4,5,-trisphosphate(IP3) ELISA kit
Reactivity	Rat
Range	5 pg/mL-1000 pg/mL
Target	IP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue_x000D_ homogenates, cell culture supernates, cell lysates
Sensitivity	2.5 pg/mL
Expiration Date	Please enquire.