

Product Datasheet

Rat SRD5A2 ELISA Kit (orb410572)

Description	Rat 3-oxo-5-alpha-steroid 4-dehydrogenase 2(SRD5A2) ELISA kit
Reactivity	Rat
Range	6 pg/mL-1200 pg/mL
Target	SRD5A2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P31214
Sensitivity	3 pg/mL
Expiration Date	Please enquire.