

Product Datasheet

Rat 8-OHdG ELISA Kit (orb410579)

Description	Rat 8-Hydroxy-desoxyguanosine,8-OHdG ELISA Kit
Reactivity	Rat
Range	0.312 ng/mL-20 ng/mL
Target	8-OHdG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.