

## Product Datasheet

### Rat PCI ELISA Kit (orb410583)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Rat typeI procollagen,PC I ELISA Kit   |
| <b>Reactivity</b>        | Rat                                    |
| <b>Range</b>             | 62.5 pg/mL-4000pg/mL                   |
| <b>Target</b>            | PCI                                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm           |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates |
| <b>Sensitivity</b>       | 15.6 pg/mL                             |
| <b>Expiration Date</b>   | Please enquire.                        |