

## Product Datasheet

### Rat Col III ELISA Kit (orb410589)

|                          |                                          |
|--------------------------|------------------------------------------|
| <b>Description</b>       | Rat Collagen Type III, Col III ELISA KIT |
| <b>Reactivity</b>        | Rat                                      |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                     |
| <b>Target</b>            | Col III                                  |
| <b>Note</b>              | For research use only                    |
| <b>Application notes</b> | Detection Wavelength: 450 nm             |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates        |
| <b>Sensitivity</b>       | 0.039 ng/mL                              |
| <b>Expiration Date</b>   | Please enquire.                          |